



ARUN
TECHNOLOGY

AURUM

THE ULTIMATE PRECIOUS
METALS ANALYSER

ACCESSORIES



SMALL SAMPLE ADAPTOR

A soft tipped clip that can be flexibly positioned to optimise the positioning of small, delicate jewellery items inside the sample chamber. Magnetic base for securing within the sample chamber.



LARGE SAMPLE ADAPTOR

A soft tipped clip mounted on a flexible positioning adaptor used to optimise the positioning of larger jewellery items and hold them securely in place for the best analysis. Magnetic base for secure mounting in the sample chamber.

